

Technical Information Sheet Product Number 9964 3M[™] Polyester Diagnostic Tape

Effective: February 2005 Supersedes: May 14, 2004

Product Description:

3MTM #9964 Polyester Diagnostic Tape consists of a clear polyester film coated on one side with a pressure sensitive acrylate adhesive and supplied on a clear silicone coated polyester release liner. Tape is provided on 6" cardboard cores with the release liner side "IN".

Composition:

Tape Layer	Description	Caliper
Backing	Clear polyester film	2.0 mils (0.051 mm) typical
Adhesive	Pressure sensitive acrylate	0.5 mils (0.013 mm) typical
Release Liner	Silicone coated polyester film	2.0 mils (0.051 mm) typical

Product Properties¹:

Adhesion to Steel ² (Reference: ASTM D3330M, 180° Adhesion Test)	Minimum 15 oz/in (425 g/25 mm)
Liner Release ² (Reference: 3M TS-239 or TS-2059, 180° Liner Release Test)	Maximum 20 g/in (20 g/25 mm)
Tape Caliper (without release liner) ² (Reference: 3M TS-303, Caliper of Materials)	2.5 ± 0.3 mils $(0.06 \pm 0.008$ mm)

Properties listed based on limited data and subject to change.

Product Dimensions:

Roll Width	Maximum 48 in (122 cm) ± 1/16 in (1.6 mm)	
Roll Length	Maximum 350 yd (320 m) ± 2%	
Splice Frequency	Maximum 2 per 350 yard (320 m) roll. (Intended for continuous processing)	
NOTE:	Product supplied on 6-inch (15.2 cm) diameter cardboard cores with liner wound "IN".	

Recommended Storage Conditions & Shelf Life:

Product as supplied in original packaging will maintain stated test properties for a period of two years from date stamped on shipping container when stored at a temperature between 50-80°F (10-27°C) and a relative humidity between 40-60 percent.

Please see reverse side for important product, safety and warranty information.

² Lot specific results may be included on a Certificate of Analysis upon request at time of order.

Product and Safety Information: User is solely responsible for determining the suitability of 3M samples and products for the intended use including any necessary safety or toxicity assessment. 3M will provide Material Safety Data Sheets and results of toxicity testing upon request. In every case before using any product in full scale production users should conduct their own tests to determine to their own satisfaction whether the product is of acceptable quality and is suitable for their particular purposes under their own operating conditions.

Notice: Nothing contained herein shall be construed to imply the nonexistence of any relevant patents or to constitute a permission, inducement or recommendations to practice any invention covered by any patent, without authority from the owners of this patent.

Warranty Information

All statements, technical information and recommendations herein are based on tests 3M believes to be reliable, but the accuracy or completeness thereof is not guaranteed. 3M warrants only that products will meet 3M's specifications at the time of shipment to the customer. 3M does not offer any other warranty and does not warrant the performance, safety or such other characteristics of Products in combination with other materials. 3M specifically DOES NOT warrant Products for any intended or unintended uses (whether or not foreseeable); for compatibility or suitability with other components or compatibility with any methods of manufacture or conversion. The foregoing warranty is made in lieu of all other warranties, expressed or implied, including the implied warranties of merchantability, fitness for a particular purpose and freedom from non-infringement.

Limitation of Remedies: If products are proven not to meet 3M's specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Products which are proven out of specification or to refund the purchase price paid for Products.

Limitations of Liabilities: The remedies provided herein are exclusive remedies against 3M for any alleged or actual nonconformance to specifications or defect or other failure in products or for 3M's performance of its supply obligations. Under no circumstances is 3M liable for any direct, incidental, special or consequential damages (including lost profits) in any way related to the product under any theory of law including, but not limited to, negligence and strict liability.

Ordering Information

Call us at Customer Service to place an order: 800-742-1994 (U.S.).

To have a sales representative contact you, to request samples or clinical and safety summaries, please contact us at 3M HELPLINE 800-228-3957 (U.S.) or 651-604-5599 (International).

For international inquiries, please contact your local country representative; our 3M Medical Specialties subsidiary contacts are listed below.

Visit our web site: www.3M.com/medicalspecialties for product and services information, news about conferences we will be attending, new product highlights or to make a direct inquiry.

International Locations

Argentina

3M Argentina S.A.C.I.F.I.A. Tel.: 54-11-4469-8200 Fax: 54-11-4469-8318

Asia Pacific

3M Asia Pacific Pte. Ltd. Singapore Tel: 65-454-8611 Fax: 65-458-5432

Australia

3M Health Care Tel: 61-2-9875-6333 Fax: 61-2-9875-6416

Belgium 3M Belgium N.V./S.A.

Hermeslaan 7 B-1831 Diegem, Belgium Tel: 32-2-722-5111 Fax: 32-2-720-0225

Brazil

3M DO Brasil Ltda. Tel: 55-19-864-7171 Fax: 55-19-864-7848

Czech Republic

3M Cesko s.r.o. Vyskocilova 1 140 00 Praha 4 Czech Republic Tel: +420-261 380 111

Fax: +420-261 380 110

Canada

3M Canada Post Office Box 5757 London, Ontario N6A 4T1 Canada Tel.: 800 364-3577 Fax: 519-452-9065

China

3M China Ltd. Tel: 86-21-6275-3535 Fax: 86-21-6275-2343

Denmark

3M a/s Fabriksparken 15 2600 Glostrup, Denmark Tel: 45-43-48-0100 Fax: 45-43-96-8596

Laboratories 3M Santé Tel: 33-1-3031-8376 Fax: 33-1-3031-8378

Germany

3M Medica Hammfelddamm 11 41453 Neuss, Germany Tel: 49 2131 14 40 00 Fax: 49 2131 14 49 99

Hungary

3M Hungária Kft. H-1138 Budapest. Váci út 178 Hungary Tel: (36-1) 270-7777 Fax: (36-1) 320-0951

India

3M India Ltd. Bangalore Raheja Paramount, 13, Residency Road, Bangalore 560 025 India Tel: 91-80-2231414 Fax: 91-80-2231450

Italy

3M Italia SPA Via San Bovio 3 -Loc. San Felice 20090 Segrate, Italy Tel: 39-02 7035-2406 Fax: 39-02 7035-3210

Japan

3M Health Care Ltd. Tel: 81-3-3709-8111 Fax: 81-3-3709-8754

Kazakhstan

3M Representation Office 153 Abay Avenue Suite 31/32 Almaty 480009 Kazakhstan Tel: +7 3272 509 944 Fax: +7 3272 509 573

Mexico

3M Mexico, S.A. de C.V. Tel: 52-55-52-70-04-00 Fax: 52-55-52-70-22-99

Middle East

3M Gulf Ltd. Dubai - UAE Tel: 971-4-265-2121 Fax: 971-4-262-7478

Netherlands

3M Nederland B.V. Tel: 31-715-450-450 Fax: 31-715-450-212

Poland

3M Poland Sp.z o.o. Kajetany k/Warszawy 05-830 Nadarzyn, Poland Tel: (48 22) 739 60 00 Fax: (48 22) 739 60 01

Portugal

3M Portugal Rua Conde de Redondo, 98 1169-009 Lisboa, Portugal Tel: 351-21-3134500 Fax: 351-21-3134680

Russia

3M Russia 125445 Moscow, Russia Ulitsa Smolnaya 24/D Meridian Commercial Tower Tel: 7 (095) 784-74-74 Fax: 7 (095) 784-74-75

Spain

3M España, S.A. Tel: 43-1-321-60-00 Fax: 43-1-321-60-02

Sweden

3M Svenska AB Bollstanäsvägen 3 SE-191 89 Sollentuna, Sweden Tel: +46-8-92 21 00

Fax: +46-8-92 24 22

3M Medical Specialties

3M Center, Building 275-5W-05 St. Paul, MN 55144-1000 USA Tel.: 800 228-3957 www.3M.com/medicalspecialties

Switzerland & Eastern

Europe 3M (Schweiz) AG Eggstrasse 93

8803 Rüschlikon, Switzerland

Tel: +41-1 724 92 44 Fax: +41-1 724 94 80

Turkev

3M Turkey Nispetiye Cad. Akmerkez Blok 3 K.5 Etiler 34340 Istanbul, Turkey Tel: 9 (0212) 3507777 Fax: 9 (0212) 2821742

3M Representation Office 77 V. Vasylkivska St., 03150 Kiev, Ukraine Tel: +380 44 490 57 77 Fax: + 380 44 490 57 75

United Kingdom & Ireland

3M Health Care Ltd. Tel: 44-1509-611-611 Fax: 44-1509-237-288

