

# Technical Information Sheet Product Number 1503 3M™ Tan Polyethylene Single Coated Medical Tape

Effective: June 2013

Supersedes all previous specifications

### **Features and Benefits:**

- Hypoallergenic Pressure Sensitive Adhesive

- Conformable

- Easy to Apply

- Provides Barrier from External Contamination

- ETO & Gamma Compatible
- Good Adhesion to Skin
- Opaque
- Warm Tan Color

## **Composition:**

Tape Caliper	5 mil (0.13 mm) Tan Polyethylene Film Tape
Backing	3.6 mil (0.09 mm) Tan Polyethylene
Adhesive	Tackified Acrylic, designed for medical/surgical applications
Release Liner	63 lb. Poly-coated Kraft paper, silicone on one side (4.9 mils/0.12 mm)

## **Typical Properties**<sup>+</sup>:

Adhesion to Steel, 180° Peel	41 oz/inch width (11.3 N/25 mm width) (1.2 Kg/25 mm width)
Liner Release	20 gm/inch width (20 gm/25 mm width)

<sup>+3</sup>M test methods are available upon request.

## **Roll Description:** Tape supplied on 3 inch (76 mm) diameter cores.

Length, Maximum	500 yards (457 m) ± 2%
Width, Maximum	48 inch (122 cm) $\pm$ 1/32 in (1.0 mm)

**Splices:** Splices on the tape are intended for continuous processing.

Roll Width	Maximum Number of Splices per Roll
≤ 8 inch (20 cm)	2 per 100 yd (91 m)
	3 per 200 yd (183 m)
	8 per 500 yd (457 m)

Information on the splice frequency in roll widths greater than 8 inch (20 cm) is available upon request.

## **Packaging and Recommended Storage:**

Product as supplied in original packaging will maintain stated properties for a period of two years from date stamped on shipping container when stored at temperatures of 50-80°F (10-27°C) and a relative humidity between 40-60 percent.

Please see reverse side for important product, safety and warranty information.

**Product and Safety Information:** User is solely responsible for determining the suitability of 3M samples and products for the intended use including any necessary safety or toxicity assessment. 3M will provide Material Safety Data Sheets and results of toxicity testing upon request. In every case before using any product in full scale production users should conduct their own tests to determine to their own satisfaction whether the product is of acceptable quality and is suitable for their particular purposes under their own operating conditions.

**Notice:** Nothing contained herein shall be construed to imply the nonexistence of any relevant patents or to constitute a permission, inducement or recommendations to practice any invention covered by any patent, without authority from the owners of this patent.

#### **Warranty Information**

All statements, technical information and recommendations herein are based on tests 3M believes to be reliable, but the accuracy or completeness thereof is not guaranteed. 3M warrants only that products will meet 3M's specifications at the time of shipment to the customer. 3M does not offer any other warranty and does not warrant the performance, safety or such other characteristics of Products in combination with other materials. 3M specifically DOES NOT warrant Products for any intended or unintended uses (whether or not foreseeable); for compatibility or suitability with other components or compatibility with any methods of manufacture or conversion. The foregoing warranty is made in lieu of all other warranties, expressed or implied, including the implied warranties of merchantability, fitness for a particular purpose and freedom from non-infringement.

**Limitation of Remedies:** If products are proven not to meet 3M's specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Products which are proven out of specification or to refund the purchase price paid for Products.

Limitations of Liabilities: The remedies provided herein are exclusive remedies against 3M for any alleged or actual nonconformance to specifications or defect or other failure in products or for 3M's performance of its supply obligations. Under no circumstances is 3M liable for any direct, indirect, incidental, special or consequential damages (including lost profits) in any way related to the product under any theory of law including, but not limited to, negligence and strict liability.

## **Ordering Information**

Call us at Customer Service to place an order: 800-742-1994 (U.S.). Visit our website: **www.3M.com/medicalspecialties** for product and services information, news about conferences we will be attending, new product highlights or to make a direct inquiry.

To have a sales representative contact you, to request samples or clinical and safety summaries, please contact us at 3M HELPLINE 800-228-3957 (U.S.) or for international inquiries, please contact your local country representative. Our 3M Medical Specialties subsidiary contacts are listed below.

## **International Locations**

Argentina 3M Argentina S.A.C.I.F.I.A. Tel.: 54-11-4339-2400 Fax: 54-11-4339-2640
Australia 3M Health Care Tel: 1 300 363 878 Fax: 1 800 060 888
Belgium

Tel: 32-2-722-5111 Fax: 32-2-720-0225	
Brazil 3M do Brasil I tda	

Tel: 55-19-3838-6661

Fax: 55-19-3838-7449

3M Belgium N.V./S.A.

Canada	
3M Canada	
Tel: 800 364-3577	

Fax: 800 341-4630

**Chile** 3M Chile S.A. Tel: (56-2) 4103000

Fax: (56-2) 4103400

## **China**3M China Ltd. Tel: 86-21-6275-3535 Fax: 86-21-5208-2205

**Colombia** 3M Colombia S.A. Tel: 57-1 4101666 Fax: 57-1 4161677

Czech Republic 3M Cesko s.r.o. Tel: +420-261 380 111 Fax: +420-261 380 110

3M a/s Tel: 45-43-48-0100 Fax: 45-43-96-8596

Denmark

## Dominican Republic

3M Dominicana Tel: (809) 530-6560 X499 Fax: (809) 530-2960

#### Finland Suomen 3M Oy Tel.: 358-9-52-521 Fax: 358-9-512-29-44

**France**Laboratoires 3M Santé
Tel: 33-1-3031-8376
Fax: 33-1-3031-8378

#### Germany 3M Medica Tel: 49 2131 14 40 00 Fax: 49 2131 14 49 99

Hong Kong 3M Hong Kong Ltd. Tel: (852) 2806 6111 Fax: (852) 2234 6044

#### **Hungary** 3M Hungária Kft. Tel: (36-1) 270-7777 Fax: (36-1) 320-0951

India 3M India Ltd. Tel: +91-80-2223 1414 Fax: +91-80-22231431

#### Indonesia PT 3M Indonesia Tel: (62-21) 520 3401 Fax: (62-21) 520 3106

Israel 3M Israel Ltd. Tel: 972-9-9615000 Fax: 972-9-9615050

3M Italia SPA Tel: +39 2 7035-2406 Fax: +39 2 7035-2484 **Japan** 3M Health Care Ltd. Tel: 81-3-3709-8289 Fax: 81-3-3709-8754

#### **Kazakhstan** 3M Representation Office Tel: +7 3272 509 944 Fax: +7 3272 509 573

**Korea** 3M Korea Ltd. Tel: 82-2-3771-4286 Fax: 82-2-786-2825

## Malaysia Sdn. Bhd. Tel: 6-03-7806 2888 Fax: 6-03-7806 2902

**Mexico** 3M México, S.A. de C.V. Tel: +52-55-5270-0400 Fax: +52-55-5270-0433

#### **Middle East** 3M Gulf Ltd. Tel: 971-4-3670-777 Fax: 971-4-3670-699

The Netherlands 3M Nederland B.V. Tel: 31-715-450-450 Fax: 31-715-450-212

#### New Zealand 3M New Zealand Tel: 64-9-477 4040

Fax: 64-9-477 6692

**Peru** 3M Peru S.A. Tel: 511-224-2728 Fax: 511-224-3171

**Philippines** 3M Philippines Tel: (632) 8133781 to 95 Fax: (632) 8145873 to 74 Poland 3M Poland Sp.z o.o. Tel: (48 22) 739 60 00 Fax: (48 22) 739 60 01

**Portugal** 3M Portugal Tel: 351-21-3134500 Fax: 351-21-3134680

Russia 3M Russia Tel: 7 (495) 784-74-74 Fax: 7 (495) 784-74-75

Singapore 3M Technologies (S) Pte. Ltd. Tel: 65-6450-8866 Fax: 65-6455-2130

**Sri Lanka** 3M Lanka (Pvt.) Ltd. Tel: (94) 11 2785701 Fax: (94) 11 4410083-4

**Spain** 3M España, S.A. Tel: 34-91-321-6000 Fax: 34-91-321-6002

Sweden 3M Svenska AB Tel: +46-8-92 21 00 Fax: +46-8-92 24 22

#### Switzerland & Eastern Europe 3M (Schweiz) AG

Tel: +41 44 724 92 44 Fax: +41 44 724 94 80

## Thailand

3M Thailand Ltd. Tel: (66)2 260 8577 Ext.175 Fax: (66)2 261 7535

Taiwan
3M Taiwan Ltd.
Tel: 886-2-2704-9011

Fax: 886-2-2706-0355 **Turkey**3M Turkey

Tel: 9 (0212) 3507777 Fax: 9 (0212) 2821742

Ukraine 3M Ukraine Tel: +380 44 490 57 77 Fax: +380 44 490 57 75

United Kingdom & Ireland 3M Health Care Ltd. Tel: 44-1509-611-611 Fax: 44-1509-237-288

#### Venezuela

3M Manufacturera Venezuela Tel: 58-241-8391911 Fax: 58-241-8391928

Viet Nam 3M Viet Nam Ltd. Tel: 84 8 640986 Fax: 84 8 641024



**3M** 

3M Medical Specialties

3M Center, Building 275-5W-05 St. Paul, MN 55144-1000 USA Tel.: 800 228-3957

FM 1455