Medical Specialties 3M Health Care 3M Center St. Paul, MN 55144-1000 651 733 1110



## PRODUCT SAFETY SUMMARY NO. 7331 (formerly MSX 6027) Gamma Stable Loop with Adhesive Effective: October 2005

## NO. 7331 GAMMA STABLE LOOP with ADHESIVE

This is a gamma-stable loop material provided in roll form. It is intended to be used on its own or in conjunction with # 7333 (MSX 6025), gamma-stable hook with adhesive. **This product is not intended for direct skin contact.** 

The adhesive used on this product has been subjected to the following biocompatibility studies.

## In Vitro Cytotoxicity (Agar Overlay)

Protocol reference: Guess, W. L. et al; "Agar Diffusion Method for Toxicity Screening of Plastics on Cultured Cell Monolayers" J. Pharm. Sci. 54:1545-1547 (1965).(HRFE 00915) **Results**: 0.0/0.0

## Acute Primary Skin Irritation in Albino Rabbits

Protocol reference: Draize: Appraisal of the Safety of Chemicals in Food, Drugs and Cosmetics (1965). Published by the Editorial Committee of the Association of Food and Drug Officials of the United States. (HRFE 00915.) **Results**: 0.0/8.0

When used as per recommended application, this product(s) is not deemed toxic, corrosive, or an irritant.

It is the responsibility of our customers to determine the final suitability of our products for their application.